

CM 23900 – TELLURITE CEFIXIME SUPPLEMENT

INTENDED USE

For isolation of *Escherichia coli* 0157:H7.

COMPOSITION

Ingredients	Concentration
Potassium tellurite	1.250mg
Cefixime	0.025mg

(per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix gently to dissolve the contents completely. Aseptically add the contents to 495 ml of sterile, molten, cooled (45-50°C) MacCONKEY SORBITOL AGAR BASE / MacCONKEY SORBITOL AGAR BASE (ISO 16654-2001). Mix well and pour into sterile petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative Mallier GmbH Buckhorn 10 41113 Aachen, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only