

CM 23902 – TOLUIDINE BLUE (0.1 g/vl)

INTENDED USE

For detection of deoxyribonuclease activity.

COMPOSITION

Ingredients	Concentration
Toluidine blue	0.1g

(per vial sufficient for 1000 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial in 5 ml of distilled water and add in 1000 ml of rehydrated DNASE TEST AGAR BASE (W/O INDICATOR). Heat to boiling to dissolve the medium completely. Sterilize by autoclaving the medium at 15 psi pressure (121°C) for 15 minutes. Cool and pour in sterile petri plates. Alternatively, 0.1% aqueous solution of toluidine blue can be used to flood the growth on DNase Test Agar plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedWet GmbH Seckestraße 11C 48153 Münster, Germany	 European Conformity	 QR Code	 Consults instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only