

CM 23909 – UREA 40% (5 ml/vl)

INTENDED USE

For detection of urease activity.

COMPOSITION

Ingredients	Concentration
Urea	2g
Distilled water	5ml

INSTRUCTION FOR USE (per vial sufficient for 100 ml medium)

Warm up the refrigerated Urea Solution to room temperature and aseptically add 5 ml in 95 ml sterile, molten, cooled (45-50°C) DIAGNOSTIC STUART'S UREA BROTH BASE (UREA BROTH BASE) / DIAGNOSTIC STUART'S UREA BROTH BASE (UREA BROTH BASE) (VEG.) / UREA BROTH BASE (CHRISTENSEN) / UREA AGAR BASE (CHRISTENSEN) (AUTOCLAVABLE) / UREA AGAR BASE (CHRISTENSEN) (AUTOCLAVABLE) (VEG.) / UREA AGAR BASE (CHRISTENSEN) / UREA AGAR BASE (CHRISTENSEN) (IS:5887 (Part-I)-1976) / MIU MEDIUM BASE / MIU MEDIUM BASE (VEG.) / HEMMES MEDIUM or 25ml in 975 ml KOHN TWO TUBE MEDIUM NO. 1 / YERSINIA IDENTIFICATION BROTH BASE. Mix well and dispense in sterile tubes.

STORAGE

Vial should be stored in sealed container at 2-8°C.

