

CM 23915 – VANCOMYCIN SUPPLEMENT

INTENDED USE

For isolation of Cronobacter sakazakii.

COMPOSITION

Ingredients	Concentration
Vancomycin	3mg

(per vial sufficient for 500 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial with 5 ml of sterile distilled water and aseptically add 2 vials to 1000 ml of sterile molten CRONOBACTER SELECTIVE BROTH (CSB) (ISO 22964:2017) / VRE AGAR BASE (VANCOMYCIN RESISTANT ENTEROCOCCI AGAR) along with 1 vial of Meropenem Supplement to Mix well and pour into sterile Petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative Mallin+Gottlieb BioScience AG 48147 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only