

CM 23918 – VITAMINS AMINO GROWTH SUPPLEMENT (Vitamins & Amino Acids Mixture)

INTENDED USE

For cultivation of wide variety of microorganisms.

COMPOSITION

Ingredients	Concentration
Vitamin B12	0.100mg
L-Glutamine	100mg
Adenine sulphate	10mg
Guanine hydrochloride	0.300mg
p-Aminobenzoic acid (PABA)	0.130mg
L-Cystine	11mg
NAD (Coenzyme I)	2.500mg
Coccarboxylase	1mg
Ferric nitrate	0.200mg
Thiamine hydrochloride	0.030mg
Cysteine hydrochloride	259mg
Dextrose	1mg
Distilled water	10ml

(per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Aseptically add this to 240 ml of sterile, molten, cooled (45-50°C) GC AGAR BASE / GC AGAR BASE (VEG.) / THAYER MARTIN MEDIUM BASE / TM 064 - CHOCOLATE AGAR BASE along with 250 ml of sterile 2% haemoglobin solution or to 500 ml of sterile, molten, cooled (45-50°C) TRANSGROW MEDIUM BASE / MYCOPLASMA UROGENITAL MEDIUM BASE or to 1000 ml sterile, molten, cooled (45-50°C) TELLURITE BLOOD AGAR BASE. Mix gently and pour into sterile petri plates.



STORAGE

Vial should be stored in sealed container at 2-8°C.

GMP Good Manufacturing Practices Certified	IVD For In Vitro Diagnostic Use	QTY. Quantity	LOT/ B. NO. Lot / Batch Number	REF Catalogue Number	 Manufacturer
 Temperature Unit	EC REP Authorized Representative Mallatier GmbH Buckhorn 10, 40142 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only