

CM 23920 – WILLIS AND HOBBS SUPPLEMENT

INTENDED USE

For isolation of Clostridium perfringens and Clostridium botulinum species from food.

COMPOSITION

Ingredients	Concentration
Neomycin sulphate (per vial sufficient for 475 ml medium)	33mg

INSTRUCTION FOR USE

Rehydrate the contents of one vial aseptically with 5.0 ml of sterile distilled water and aseptically add it to 400 ml of sterile, molten, cooled (45-50°C) WILLIS AND HOBBS MEDIUM BASE (IS: 5887 (Part 4) : 1999) along with 15 ml of Egg Yolk Emulsion and 60 ml of sterile skim milk.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedWet GmbH Sankt Augustin 10 50153 Aachen, Germany	 CE European Conformity	 QR Code	 Consults instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only