

CM 23922 – YEAST AUTOLYSATE SUPPLEMENT

INTENDED USE

For the cultivation of *Neisseria* species.

COMPOSITION

Ingredients	Concentration
Yeast autolysate	5g
Dextrose	0.500g
Sodium bicarbonate	0.075g

(per vial sufficient for 500 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 10 ml sterile distilled water. Mix well & aseptically add to sterile, molten, cooled (45-50°C) 490 ml of GC AGAR BASE / GC AGAR BASE (VEG.) containing Haemoglobin along with GC supplement w/ antibiotics or in CHOCOLATE AGAR BASE along with 5% v/v defibrinated blood (after addition of blood mix and heat at 80°C for 10 mins). Mix well and pour into sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MedWet GmbH Seestraße 11 40113 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only