

CM 23924 – β LACTAMASE I SUPPLEMENT (W/ 3300 IU OF PENICILLINASE)

INTENDED USE

For inactivation of the penicillin group of antibiotics.

COMPOSITION

Ingredients	Concentration
Penicillinase	3300 IU of activity

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial with 5 ml of distilled water. After reconstitution, filter sterilize the solution with 0.2 micron syringe filter. Add appropriate amount of solution depending on the application. Remaining sterilized solution can be stored at 2-8° C for 4 weeks.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedNet GmbH Barkstrasse 10, 48163 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only

