

## CM 23707 – ALBUMIN GLUCOSE SUPPLEMENT

### INTENDED USE

For rapid cultivation of Mycobacterium tuberculosis.

### COMPOSITION

Ingredients	Concentration
Albumin fraction V	1g
Dextrose	1.500g
Distilled water	20ml

(Per vial sufficient for 500ml medium)

### APPLIED IN MEDIUM


DUBOS BROTH BASE

### INSTRUCTION FOR USE

Warm up the refrigerated supplement to room temperature. Aseptically add to sterile, cooled (45-50°C) 180 ml Dubos Broth Base. Mix well and pour into sterile screw capped tubes.

### STORAGE:

Vials should be stored in sealed container at 2-8°C.



NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

\*For Lab Use Only

