

CM 23744 – CHLORTETRACYCLINE SELECTIVE SUPPLEMENT

INTENDED USE

For selective isolation of *Penicillium viridicatum* and related species from foods.

COMPOSITION

Ingredients	Concentration
Chlortetracycline	50mg

(Per vial sufficient for 1000 ml medium)

INSTRUCTION FOR USE

Rehydrate the content of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add it to 995 ml of sterile, molten, cooled (45-50°C) APRY BROTH BASE / PENTACHLORO ROSE BENGAL YEAST EXTRACT AGAR BASE (PRYES AGAR) (as per APHA) / ROSE BENGAL AGAR W/ CHLORTETRACYCLINE / SABOURAUD GLUCOSE AGAR W/ANTIBIOTICS / SABOURAUD DEXTROSE AGAR MEDIUM W/ ANTIBIOTICS (as per IP) / AGAR MEDIUM C (SABOURAUD-GLUCOSE AGAR WITH ANTIBIOTICS) (as per EP/BP) / YEAST MOULD CHLORAMPHENICOL AGAR, MODIFIED. Mix well and dispense as desired.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MedWet GmbH Seckestraße 15, 40162 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only

