

CM 23769 – DOYLE'S ANTIBIOTIC SUPPLEMENT

INTENDED USE

For selective isolation of Campylobacter species.

COMPOSITION

Ingredients	Concentration
Vancomycin	7.500mg
Trimethoprim	2.500mg
Polymyxin B sulphate	10000IU
Cycloheximide	25mg

(per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 5 ml of 50% ethanol. Mix well and aseptically add to 460 ml of sterile, molten, cooled (45-50°C) DOYLES ENRICHMENT BROTH BASE along with 35 ml sterile lysed horse blood. Mix well and dispense as desired.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative Machner GmbH Buckhorn 10, 40143 Essen, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only