

CM 23778 – FRASER ENRICHMENT SUPPLEMENT

INTENDED USE

For isolation and enrichment of *Listeria monocytogenes* from food and environmental samples.

COMPOSITION

Ingredients	Concentration
Acriflavine hydrochloride	25mg
Nalidixic acid	4mg

(per vial sufficient for 1000 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of one vial aseptically with 10 ml of sterile distilled water. Mix well and aseptically add it to 990 ml sterile, cooled (45-50°C) FRASER SECONDARY ENRICHMENT BROTH BASE. Mix well and dispense as desired.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 Best Before	 QTY. Quantity	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 LOT/ B. NO. Lot / Batch Number	 Consults Instructions for Use	 QR Code	

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only