

CM 23789 – GROWTH SUPPLEMENT I FOR MSM

INTENDED USE

A growth supplement for the enrichment of Salmonella species.

COMPOSITION

Ingredients	Concentration
Vitamin B12	0.100mg
Distilled Water	1ml


(per vial sufficient for 1000 ml medium)

INSTRUCTION FOR USE

Warm up the refrigerated contents of one vial to room temperature. Aseptically add the contents of one vial of Growth Supplement I for MSM (TS 299) to 1000 ml sterile, cooled (45-50°C) MSM BROTH BASE along with Growth Supplement II for MSM (TS 300). Mix well and dispense as desired.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedWet GmbH Buckhorn 10 40162 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only