

CM 23809 – LANDER VPT SUPPLEMENT A

INTENDED USE

For selective enrichment and transport of *Campylobacter* spp.

COMPOSITION

Ingredients	Concentration
Polymyxin B	10000IU/L
Trimethoprim	20mg/L
Vancomycin	40mg/L


(Per vial sufficient for 1000ml medium)

INSTRUCTION FOR USE

Rehydrate the contents aseptically of 1 vial into to 1000 ml of sterile, molten LANDER TRANSPORT & STORAGE MEDIUM. Mix carefully and dispense as desired.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative Mallatier GmbH Buckhorn 10, 40112 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only