

CM 23814 – LEGIONELLA GROWTH SUPPLEMENT (BCYE)

INTENDED USE

Forenhancinggrowth of Legionella species.

COMPOSITION

Ingredients	Concentration
ACES buffer/ Potassium hydroxide	5g
Ferric pyrophosphate, soluble	0.125g
L-cysteine hydrochloride	0.20g
alpha-Ketoglutarate	0.50g
Distilled water	50ml

(Per vial sufficient for 500 ml medium)

INSTRUCTION FOR USE

Warm up the refrigerated supplement to 45-50°C, shake well to form a uniform suspension. Avoid frothing. Aseptically add to 440 ml sterile, molten, cooled (45-50°C) LEGIONELLA AGAR BASE / LEGIONELLA AGAR BASE W/O CHARCOAL. If desired, aseptically add the rehydrated contents of 1 vial of Legionella (GVPC) Selective Supplement. Mix well and pour into sterile petri plates. The final pH of the medium will be 6.9 ± 0.2. In case of non-incorporation of Legionella (GVPC) Selective Supplement, aseptically add 10 ml sterile distilled water to bring the total volume to 500 ml medium.

Note: On storage at 2 to 8°C; due to the combination of ingredients, slight fine precipitate may develop which may be evenly distributed before addition to medium.

STORAGE

Vialsshouldbe stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MedWet GmbH Buckhorn 15, 40142 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For LabUse Only

