

CM 23817 – LEGIONELLA SELECTIVE SUPPLEMENT (GVPN)

INTENDED USE

For enhancing growth of Legionella species.

COMPOSITION

Ingredients	Concentration
Glycine	1.500g
Vancomycin hydrochloride	0.500mg
Polymyxin B sulphate	40000IU
Natamycin	20mg

(Per vial sufficient for 500 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 10 ml of sterile distilled water. Mix well and aseptically add to 440 ml of sterile, molten, cooled (45-50°C) LEGIONELLA AGAR BASE / LEGIONELLA AGAR BASE W/O CHARCOAL along with 50 ml of Legionella Growth Supplement BCYE. Mix well and pour into sterile petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedWet GmbH Buckhorn 16, 40143 Essen, Germany	 CE European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For LabUse Only