

CM 23821 – LEGIONELLA SELECTIVE SUPPLEMENT-IV (MWY)

INTENDED USE

For selective isolation of Legionella species.

COMPOSITION

Ingredients	Concentration
Polymyxin B sulphate	25000Unit
Glycine	1.500g
Anisomycin	40mg
Vancomycin	0.500mg
Bromo thymol blue	5mg
Bromo cresol purple	5mg

(Per vial sufficient for 500 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 5 ml of 50% ethanol. Mix well and aseptically add to 500 ml of sterile, molten, cooled (45-50°C) BUFFERED CHARCOAL YEAST EXTRACT AGAR BASE / BUFFERED CHARCOAL YEAST EXTRACT AGAR MEDIUM (ISO 11731-2:2017) along with the rehydrated contents of 1 vial of Legionella Supplement. Mix well and pour into sterile Petri plates.

STORAGE

Vialsshouldbe stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MedWet GmbH Buckhorn 16, 40147 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For LabUse Only

