

CM 23822 – LEGIONELLA SUPPLEMENT

INTENDED USE

For enhanced growth of Legionella species.

COMPOSITION

Ingredients	Concentration
L-Cysteine hydrochloride	200mg
Ferric pyrophosphate, soluble	125mg
Distilled water	5ml

(Per vial sufficient for 500 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial with 5 ml sterile distilled water and aseptically add to 490 ml of sterile, molten, cooled (45-50°C) TM 684 - BUFFERED CHARCOAL YEAST EXTRACT AGAR BASE / BUFFERED CHARCOAL YEAST EXTRACT AGAR MEDIUM (ISO 11731-2:2017) / LEGIONELLA AGAR BASE along with Legionella Selective Supplement IV (MWY) (TS 018) as recommended. Mix well and pour into sterile Petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MediMet GmbH Buckstrasse 10, 48143 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For LabUse Only

