

CM 23823 – LEPTOSPIRA ENRICHMENT SUPPLEMENT (Double Pack)

INTENDED USE

For cultivation of Leptospira species.

COMPOSITION

Ingredients	Concentration
Calcium chloride	2mg
Magnesium chloride	2mg
Zinc sulphate	0.800mg
Copper sulphate	0.060mg
Ferrous sulphate	10mg
Vitamin B12	0.040mg
Bovine serum albumin	2g
Polysorbate 80	2.500ml
Glycerol	1ml
Distilled water	20ml

(Per vial sufficient for 180 ml medium)

INSTRUCTION FOR USE

Warm up the refrigerated Leptospira Enrichment to room temperature and aseptically add 1 vial to 180 ml of sterile, molten, cooled (45-50°C) LEPTOSPIRA MEDIUM BASE. Mix well and dispense as desired.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MediStar GmbH Buckhorn 10, 40142 Alsdorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For LabUse Only

