

CM 23841 – MKTT NOVOBIOCIN SUPPLEMENT

INTENDED USE

For selective enrichment and isolation of Salmonellae.

COMPOSITION

Ingredients	Concentration
Novobiocin	40mg


(Per vial sufficient for 1000ml medium)

INSTRUCTION FOR USE

Rehydrate contents of 1 vial aseptically with 5 ml of sterile distilled water and aseptically add to sterile, cooled (45-50°C) MULLER KAUFFMAN TETRATHIONATE NOVOBIOCIN BROTH BASE (ISO 6579-1:2017). Mix well and dispense as desired.

STORAGE

Vials should be stored in sealed container at 2-8°C.



NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.