

CM 23867 – PEMBA SUPPLEMENT

INTENDED USE

For selective cultivation of *Bacillus cereus*.

COMPOSITION

Ingredients	Concentration
Polymyxin B sulphate	10000Unit


(per vial sufficient for 90 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of one vial aseptically with 5 ml sterile distilled water. Mix well and aseptically add it to 90 ml of sterile, molten, cooled (45-50°C) POLYMYXIN PYRUVATE EGG YOLK MANNITOL BROMOTHYMOL BLUE AGAR BASE (PEMBA) along with 5 ml of Egg Yolk Emulsion. Mix well and pour into sterile petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 Best Before	 QTY. Quantity	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 LOT/ B. NO. Lot / Batch Number	 Consults Instructions for Use	 QR Code	

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only