

CM 23885 – S.T.A. SELECTIVE SUPPLEMENT

INTENDED USE

Recommended for the selective isolation of *Brochothrix thermosphacta* from meat products.

COMPOSITION

Ingredients	Concentration
Streptomycin sulphate	250mg
Cycloheximide	25mg
Thallous acetate	25mg


(per vial sufficient for 500 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial with 5 ml of sterile distilled water. Mix well and aseptically add it to 495 ml of sterile, molten, cooled (45-50°C) S.T.A. Agar Base. Mix well and pour into sterile petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 Best Before	 QTY. Quantity	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 LOT/ B. NO. Lot / Batch Number	 Consults Instructions for Use	 QR Code	

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only