

CM 23903 – TRIBUTYRIN (10ml/vl)

INTENDED USE

For detection of lipolytic microorganisms.

COMPOSITION

Ingredients	Concentration
Tributyrin	10ml


(per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Warm up the refrigerated solution to room temperature. Add contents of one vial aseptically to 990 ml of TRIBUTYRIN AGAR BASE W/O TRIBUTYRIN and then sterilize by autoclaving at 15 psi pressure for 15 minutes. Mix well and pour into sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 Best Before	 QTY. Quantity	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 LOT/ B. NO. Lot / Batch Number	 Consults Instructions for Use	 QR Code	

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only